

# Technical Data Sheet

Doc No– Ab/TDS/LH-73

Rev No: 03 Date 10/2021-22

## **TECHNICAL DATA SHEET**

**Abdos Labtech Pvt. Ltd.**

**Handypette Pipette Aid.**

**Valid for Catalog No. E11218, E11216, E11222**

### **1. Description/Specification**

**1.1 Description:** Handypette Pipette Aid.

1.1.1 Handypette Pipette Aid ensures quick & accurate operation.

1.1.2 Colour Coded by volume with green for 10 ml, red for 25 ml & blue for 2ml.

1.1.3 Users can easily insert glass & plastic pipette into the collar with stabilizing fingers for secure hold.

1.1.4 Resistant to acids, alkalis & easily disassembly for cleaning.

1.1.5 Rotating the knurled thumbwheel draws liquid to desired volumes. To release liquid apply light pressure to the fast release lever. For graduation dispensing rotate the thumb wheel.

### **1.2 Design:**

Handypette Pipette

### **1.3 Performance Specifications:**

<b>Code</b>	<b>Capacity (ml)</b>	<b>No. per Case</b>
<b>E11218</b>	10	4
<b>E11216</b>	25	4
<b>E11222</b>	2	4



**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India  
**An ISO 9001 ,13485, 14001 & CE Certified**  
A Division of Abdos

**Regd. Office:**  
13A, Government Place (East)  
Kolkata – 700069, India  
**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**  
**E : labtech@abdosindia.com**  
**W: www.abdoslifesciences.com**  
**CIN : U85100WB2008PTC123903**

# Technical Data Sheet

## 1.4 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release

## 1.5. Basic Features:

**1.5.1 Autoclavability:** NO, Non autoclavable

**1.5.2 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.5.3 CE certified**



**ABDOS**

**Shiv Ganga Industrial Estate,**  
Plot no.1 Lakeshwari, Bhagwanpur,  
Roorkee – 247667 Uttrakhand, India

**An ISO 9001 ,13485, 14001 & CE Certified**

A Division of Abdos

**Regd. Office:**

13A, Government Place (East)  
Kolkata – 700069, India

**GST : 05AAGCA8605A1Z7**

**P : +91 (0) 1332 235277**

**E : labtech@abdosindia.com**

**W: www.abdoslifesciences.com**

**CIN : U85100WB2008PTC123903**