

# Technical Data Sheet

Doc No Ab/TDS/CT-21

Rev No: 05 Date 02/2022-23

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Swirl Pro High Speed Micro Centrifuge**

**Valid for Catalog No. E11520**

### **1. Description Specifications:**

#### **1.1 Description: Swirl Pro High Speed Micro Centrifuge**

1.1.1 Swirl Pro is a micro centrifuge used for sedimentation of cells and virus, separation of sub- cellular organelles such as DNA, RNA, proteins & lipids.

1.1.2 Swirl Pro has a large ergonomic design with easy to read display & sound alert.

1.1.3 Both RPM & G-force can be displayed.

1.1.4 Swirl Pro adapts a Brushless DC Motor.

1.1.5 Special multi flow air-cooling design helps to keep lower rotor temperature.

1.1.6 Special multi-flow air cooling design helps to keep lower rotor temp.

1.1.7 Unique cooling system design minimizes rotor heating to less than 10°C after running for one hour

#### **1.2 Design: Swirl Pro High Speed Micro Centrifuge**

#### **1.3 Performance Specifications:**

<b>Specifications</b>	<b>E11520</b>
<b>Voltage</b>	Single-phase, AC 220-240V
<b>Frequency</b>	50-60 Hz
<b>Current</b>	5A
<b>Max Speed</b>	15000 rpm (200-15000rpm) with increment of 10 rpm
<b>Max RCF</b>	21380xg with increment of 10xg
<b>Speed Accuracy</b>	±20 rpm
<b>Rotor Capacity</b>	(i) 5ml x 18 places (ii) 0.2ml/0.5ml/1.5ml/2ml x 24 places (iii) 0.5ml x 36 places (iv) PCR 8x4 places



Abdos Labtech Private Limited  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277  
Email:labtech@abdosindia.com  
Web: [www.abdoslivesciences.com](http://www.abdoslivesciences.com)  
GST: 05AAGCA8605A1Z7

# Technical Data Sheet

<b>Run Time</b>	30sec-99min- HOLD (Continuous operation)
<b>Driving Motor</b>	Brushless DC Motor
<b>Safety Devices</b>	(i) Dual door interlock (ii) Over- speed detection (iii) Over temperature detection (iv) Automatic internal diagnosis
<b>Noise Level</b>	≤ 58 dB
<b>Additional Features</b>	(i) Speed/RCF switch (ii) Short-time run function (iii) Sound-alert function
<b>Acceleration / Breaking time</b>	25 sec for Both Acceleration & Breaking
<b>Dimension</b>	280x364x266 mm
<b>Weight</b>	12 kg
<b>Permissible ambient Temp.</b>	+4 to +40°C
<b>Permissible ambient Humidity</b>	80%

**1.4 Material:** ABS (Acrylonitrile Butadiene Styrene)

**1.5 Color:** White

**1.6 E11520** – Swirl Pro High Speed Micro Centrifuge, CE certified

**1.7 Spare Accessories**

1.7.1 **E11521** – 24 Place Rotor for 0.2, 0.5, 1.5 & 2.0ml

1.7.2 **E11522** – 18 Place Rotor for 5.0ml

1.7.3 **E11523** - 36 Place Rotor for .02 & 5.0ml

1.7.4 **E11524**- 32 (8x4) Place Rotor for 0.2 ml PCR

1.7.5 **E11525**- Rotor Adaptor for 2.0ml-0.2ml used with E11521

1.7.6 **E11526**- Rotor Adaptor for 2.0ml-0.5ml used with E11521

1.7.7 **E11527**- Rotor Adaptor for 0.5ml-0.2ml used with E11523

**Warranty:** 12 months from the date of Installation or 15 months from the date of Invoicing.



Abdos Labtech Private Limited  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakheswari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277  
Email:labtech@abdosindia.com  
Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)  
GST: 05AAGCA8605A1Z7

# Technical Data Sheet

## Instruction for Use

Storage temp. Max 30°C  
Before using the product check .  
Check all the information printed on pouches or boxes.  
Always wear gloves while handling specimen.

### WARNING

**DO NOT** store strong oxidizing agents.  
**DO NOT** place any plastic labware in a flame.  
**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure  
Consult local fire codes prior to storage of flammable liquids.



**Abdos Labtech Pvt. Ltd.**

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India  
Email : labtech@abdosindia.com, www.abdoslifesciences.com  
Ph: +(91) 11 41081897



**Abdos Labtech Private Limited**  
(A Division of Abdos Group)  
(An ISO 9001,13485,14001,18385,GMP & CE certified company)  
Shiv Ganga Industrial Estate  
Plot No-1, Lakshwari, Bhagwanpur,  
Roorkee-247667, (UK) INDIA

**Registered office**  
13A Government  
Place East, 2nd Floor Kolkata  
West Bengal – 700069  
CIN: U85100WB2008PTC123903  
**Corporate Office**  
319, 3rd Floor, DLF Tower A, DDA Dist. Centre  
Jasola, New Delhi - 110025

**Ph. No: +91- (0) 1332235277**  
**Email:labtech@abdosindia.com**  
**Web: [www.abdoslifesciences.com](http://www.abdoslifesciences.com)**  
**GST: 05AAGCA8605A1Z7**